

# Precise application in complex thin-walled anatomy

*Benefit from high-pressure needle-free injection, precise cutting and a new flexibility for retroflexion.*

**TAILOR YOUR OWN TECHNIQUE.**

**erbe**  
power your performance.

## HYBRIDknife® flex I-Type

For efficient  
3rd space endoscopy

Erbe Elektromedizin GmbH · Waldhoernlestrasse 17 · 72072 Tuebingen · Germany · +49 7071 755-0 · info@erbe-med.com

JOIN US ON     

Product not yet commercially available.





Think pressure. Think precision. Think flexibility.



### Prepare your working space with high-pressure

You can separate tissue even over folds or submucosal fibrosis. Additionally, spatial cushion provides a clear view of the operating site.

### Speed up your precise resection




Cut precise and limit the thermal damage with the short and thin electrode even in the right colon.

### Use instant coagulation for hemostasis

Several tip designs together with a variety of dissection modes facilitates a range of coagulation effects.

### Reach challenging anatomy

Even in retroflexion, the electrode remains stable and the flexibility of the probe allows it to pass smoothly through the working channel.



Marketing claim list 0188716, V005

erbe





## The right mode for every step

*HYBRIDknife® flex with the all-in-one GI workstation offers you the right mode for electrosurgical precision. endoCUT® for incision, preciseSECT for dissection and coagulation, to name just a few.*

**TAILORED FOR YOUR OWN TECHNIQUE.**

### Technical properties

	Tubing length	2300 mm
	Outer diameter (distal end)	2.6 mm
	Outer diameter (tubing)	2.1 mm
	Electrode length	1.5 mm / 2.0 mm
	Type	I
	Retro-flexibility	More than 210°
	I-Type 1.5 mm	No. 20150-112
	I-Type 2.0 mm	No. 20150-113